



Clean**up**

Antibacterial **occluder**



Ps Anesthesia[®] sa
Innovation in health

Cleanup

Antibacterial occluder

“Passive protection and disinfection device for needle-free connectors”

Without a proper disinfection, infusion ports are potential focus of intraluminal contamination.¹

1' Disinfects in 1 minute and protects for up to 7 days, provided it is not removed.



Ensures safety and effectiveness in disinfection protocol compliance.




Contributes to minimizing bloodstream infection rates.



Luer lock design, compatible with all needle-free connectors on the market ³



Presentation in die cut liner with upper orifice, facilitates the health professional an easy access to the device.



1,2.Gorski L., Hadaway L., Hagle ME., McGoldrick M., Orr M., y Daellman D. (2016). Standar of practice. Journal of Infusion Nursing; 34(1S): S73, S74
3. ISO 80369-7