

ne
nutren



Ps Anesthesia® sa
Innovation in health




“ The administration system for a safe enteral nutrition with ENFit connection ”



Connection systems constitute a key **factor in the safety of the patient.**



Inverted position of the male/female connectors in relation to the intravenous infusion systems.



Incompatibility due to design.

The difference in diameter and angles of screw closures prevents the connection of the Luer Lock system (ISO 80369-7 ex 594) used in tubes, syringes, and vascular access catheters.

**Color is used as the primary attribute for identifying function.
VIOLET = ENTERAL**

**Scan the QR code to explore
the entire catalog**

