



Comprehensive solutions
designed for prevention



Specialized nutritional support is not free of risk. Connection errors between parenteral and enteral administration systems occur and jeopardize the safety of the patient.

Nutren is the administration system for a safe enteral nutrition with ENFit connection according to the rule of design ISO 80369-3, designed by PSA+

Safe Enteral Nutrition

LATEX
FREE



SAFE CONNECTION
ENFit



DEHP
FREE

Reduces the risk of administering medication and nutritional formula via parenteral.




Connection systems constitute a **key factor in the safety of the patient.**

REDUCES THE RISK OF CONNECTION ERRORS



Incompatibility due to design: The difference in diameters and angles of screw closures prevent the connection of the Luer system (ISO 80369-7 ex594) using tubes, needles, and catheters for vascular access.



Inverted position: of the male/female connectors in relation to the intravenous infusion systems.




IMPROVES
PATIENT'S SAFETY

ISO 80369-3



A safe transition to a safe connection



Adaptador para jeringas SLIM FEMALE ENFit




Adaptador ENPlus




Adaptador MALE ENFit-ENLock

The inclusive implementation of the ENFit safe connection system involves a transition phase. **PSA+ includes ENLock adapters** which allow the connection of the existing enteral tubes and **ENPlus adapters** so as to avoid losses in the connection with the feeding tube due to the use of parenteral infusion.



Macro extension for enteral administration.

Violet PVC tubing 4.1mm, ENFit connectors, male and female




Micro extension for enteral administration.

Violet PVC tubing 2.5mm, ENFit connectors, male and female.



Enteral administration set.

Enteral administration system, ENPlus connector, drip chamber (20drops=1ml), 1 ENFit port (male), ENFit connector (female).



Enteral administration set

Multi spike Circuit for enteral administration with double entry to the drip chamber (20 drops=1ml), ENFit connector (female).



Enteral Nutrition Bottle. 125cc suitable for pasteurization.

Enteral Nutrition Bottle. 650cc suitable for pasteurization.

Enteral Nutrition Bottle. 300cc suitable for pasteurization.



DISPOSITIVOS
MEDICOS
PI 13485-2016

www.anesthesiasa.com



Administración Nacional de Medicamentos,
Alimentos y Tecnología Médica
PM 1583-18



Ps Anesthesia innovative elements in health

Anesthesia family:

Ps Anesthesia **PERISET**

Ps Anesthesia **RAQUISET**

Aneskit



ne
nutren


Cleanup

vitrolife

labkit+

hemolife

CVCmax



Pte. Perón N° 100 | Brandsen (CP 1980) | PCIA Bs. As.
Tel.: (02223) 44 4260 | Fax (02223) 44 4595